

◆◆ **PRODUCT SPECIFICATION** ◆◆



Medatro® Medical trolley B17

Patient monitor medical trolley
 Medical cart for hospital ICU emergent devices
 Professional designed hospital device

Model: B17

CONFIGURATION	PARAMETER
Model	B17
Outline size	φ 600*890mm
Column	78*100*810mm
Base	φ 600*70mm
Mounting platform size	240*225*8mm
Texture	Aluminum+steel
Color	White
Caster	4 inch*5 pcs with brake
Wheel	Silent wheels
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	12kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	15kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

