

◆◆ **PRODUCT SPECIFICATION** ◆◆



MediFocus medical trolley for
eVolution 3E-HP ventilator

CONFIGURATION	PARAMETER
Model	E01
Overall Size	580*510*920mm
Column size	135*60*730mm
Base size	580*510*50mm
Mounting platform size	487*378*73mm
Texture	Stainless steel + aluminum
Color	White + Blue
Caster	4 inch*4 pcs (brake+swivel)
Wheel	Noiseless wheels
Bearing Capacity	50kg
Max push speed	2m/s
Max tilt	15°
Net weight	22.2kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	25kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

