

◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	E05
Column	210*55*735mm
Mounting platform	352*295*133mm
Handle	Yes (250*210*5mm)
Basket	No
Humidifier hanger	No
Wire wrap	Yes (186*55*75mm)
Base	Yes (622*548*70mm)
Wheel	Yes (4X4" brake)
Bearing Capacity	50kg
Max push speed	2m/s
Max tilt	15°
Net weight	28kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	30kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

