

◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	H02
Outline	581*520*1220mm
Mounting platform	563*225*35mm
Base	581*520*50mm
Column	1000*122*35.5mm
Handle	Yes
Monitor stand	Yes
Shelf	500*490*20mm
Wheel	4X4" swivel and brake
Bearing Capacity	50Kg
Max push speed	2m/s
Max tilt	15°
Net weight	37.5kg
Carton dimension	1300*635*300mm
Gross weight	42kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

