

◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	J01
Outline size	533.4*430*710mm
Column	φ 34*617*3mm
Base	533.4*430*27mm
Mobile stand	φ 25*579*1.5mm
Wheel	4X2" swivel and brake
Bearing Capacity	5kg
Max push speed	2m/s
Max tilt	15°
Net weight	7.4kg
Carton dimension	560*450*165mm
Gross weight	8.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengen Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

