


◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	K06
Outline size	600*550*1140mm
Column	70*135*1000mm
Base	600*550*165mm
Mounting	500*430*30mm
Drawer	400*425*95mm
Monitor stand	405*35mm
Handle	Yes
Wheel	4X4" swivel and brake
Bearing Capacity	100kg
Max push speed	2m/s
Max tilt	15°
Net weight	50kg
Carton dimension	1060*570*260mm (L*W*H)
Gross weight	53kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

