


◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	K12
Column	70*135*980mm(L*W*H)
Display bracket	405*35mm
Base	600*550*165mm
Handle	Yes (ϕ 22*450*160mm)
Drawer	Yes (400*425*95mm)
Humidifier hanger	No
Wire hanger	No
Wheel	Yes (2X4" brake)
Bearing Capacity	100Kg
Max push speed	2m/s
Max tilt	15°
Net weight	49kg
Carton dimension	805*635*300mm (L*W*H)
Gross weight	57.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

