


◆ PRODUCT SPECIFICATION ◆



CONFIGURATION	PARAMETER
Model	L01
Column	200*50*935mm(L*W*H)
Mounting platform	610*500*30mm
Base	810*630*60mm
Handle	Yes (330*30*15mm)
Partition	Yes (two*550*500*30)
Drawer	Yes (470*210*30mm)
Display bracket	Yes (60*330mm)
Keyboard	Yes
Wheel	Yes (2X4" brake)
Bearing Capacity	100Kg
Max push speed	2m/s
Max tilt	15°
Net weight	55kg
Carton dimension	670*620*970mm (L*W*H)
Gross weight	70kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com