

◆◆ **PRODUCT SPECIFICATION** ◆◆



CONFIGURATION	PARAMETER
Model	M01
Column	40*710mm(ϕ *H)
Mounting platform	260*200*8mm
Handle	Yes
Basket	Yes
Base	ϕ 520*40mm
Humidifier hanger	Yes (55*40*16mm)
Wire hanger	Yes
Wheel	Yes (3X3'' universal swivel and 2x3'' brake)
Bearing Capacity	10kg
Max push speed	2m/s
Max tilt	15°
Net weight	7kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	9.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

