

Circuit hanger

## SA-03

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Column loading capacity: 5Kg

Degree: 0° -210°

Length: 1380mm

Diameter:  $\phi 14$

Plug-in:  $\phi 16*65$ mm

Thread: M14\*1.5mm

Flex tube:  $\phi 10*140$ mm

Flex tube loading capacity: 0.75Kg

Joint: 2

Weight: 1.26Kg

Carton size: 530\*95\*80mm

Gross weight: 1.44Kg

