

◆◆ PRODUCT SPECIFICATION ◆◆



MediFocus medical trolley for
Vyaire carefusion VELA ventilator

CONFIGURATION	PARAMETER
Model	B07
Outline Size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Mounting platform	282*263.5*8mm
Humidifier hanger	Yes
Handle	Yes
Oxygen bottle support	Yes
Wheel	3 inches*5 pcs (brake+swivel)
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	11.3kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	14kg