


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	B12
Column	78*100*810mm(L*W*H)
Mounting platform	230*245*32mm
Handle	Yes
Basket	Yes (330*210*135mm)
Humidifier hanger	Yes (55*40*16mm)
IV pole	Yes (Φ 19*780mm)
Wire hanger	Yes
Wheel	Yes (3X3" universal swivel and 2x3" brake)
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	10.2kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	12.8kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengen Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com