


◆◆ PRODUCT SPECIFICATION ◆◆



Medatro® Medical trolley B19

Mobile computer cart can be used in hospital

Medical trolley for hospital ICU emergency devices

Standard VESA support

Model: B19

CONFIGURATION	PARAMETER
Model	B19
Outline size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Mounting platform size	100*78*128mm
Drawer size	430*290*131.5mm
Texture	Stainless steel
Color	Medical white
Caster	4 inch*4 pcs (brake+swivel)
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	20kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	23kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tenggren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

