

◆◆ PRODUCT SPECIFICATION ◆◆



MediFocus medical trolley for

Bellavista 1000 ventilator

CONFIGURATION	PARAMETER
Model	B15
Outline Size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Mounting platform	300*230*8mm
Humidifier hanger	Yes
Handle	Yes
Oxygen bottle support	Yes
Wheel	3 inches*5 pcs (brake+swivel)
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	14.9kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	18kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com