

◆◆ PRODUCT SPECIFICATION ◆◆

MediFocus medical trolley for
Patient monitor

CONFIGURATION	PARAMETER
Model	B22
Outline Size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Mounting platform size	360*265*8mm
Texture	Aluminum+steel
Color	White
Caster	3 inch*5 pcs (brake and universal)
Wheel	Silent wheels
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	20kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	23kg