

## ◆◆ PRODUCT SPECIFICATION ◆◆



| CONFIGURATION          | PARAMETER                            |
|------------------------|--------------------------------------|
| Model                  | B24                                  |
| Outline size           | Φ 600*890mm                          |
| Column                 | 78*100*810mm                         |
| Base                   | Φ 600*70mm                           |
| Basket                 | 330*210*135mm                        |
| Humidifier hanger      | 55*40*16mm                           |
| Oxygen Cylinder Holder | Yes                                  |
| Handle                 | Yes                                  |
| Wheel                  | 3X3" universal swivel and 2x3" brake |
| Bearing Capacity       | 30kg                                 |
| Max push speed         | 2m/s                                 |
| Max tilt               | 15°                                  |
| Net weight             | 17kg                                 |
| Carton dimension       | 900*570*210mm (L*W*H)                |
| Gross weight           | 19.6kg                               |

**Beijing MediFocus Medical Co.,Ltd**

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: [info@bjmedifocus.com](mailto:info@bjmedifocus.com)

[www.bjmedifocus.com](http://www.bjmedifocus.com)

