

◆◆ PRODUCT SPECIFICATION ◆◆

CONFIGURATION	PARAMETER
Model	B26
Outline size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Basket	330*210*135mm
Humidifier hanger	55*40*16mm
Wire Hanger	Yes
Handle	Yes
Wheel	3X3" universal swivel and 2x3" brake
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	8.9kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	11.5kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tenggren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.comwww.bjmedifocus.com