


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	B28
Outline size	Φ 600*890mm
Column	78*100*810mm
Base	Φ 600*70mm
Basket	330*210*135mm
Mounting	275*200*5mm
Wire wrap	Yes
Handle	Yes
Wheel	5X4" swivel and brake
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	19.4kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	22kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

