


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	D01
Column	51*810mm (Φ *H)
Mounting platform	230*245*32mm
Handle	Yes
Basket	Yes (330*210*135mm)
Humidifier hanger	Yes (55*40*16mm)
Wire wrap	Yes (280*35*2mm)
IV pole	Yes (Φ 19*490*200mm)
Wheel	Yes (3X3" universal swivel and 2x3" brake)
Bearing Capacity	15kg
Max push speed	2m/s
Max tilt	15°
Net weight	18.5kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	21kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

 Email: info@bjmedifocus.com
www.bjmedifocus.com