


◆◆ PRODUCT SPECIFICATION ◆◆



MediFocus medical trolley for
HVJ-88oAventilatorF

CONFIGURATION	PARAMETER
Model	D04
Outline size	Φ 600*960mm
Column	Φ 51*800mm
Base	Φ 600*80mm
Mounting platform	400*314*5mm
Handle	Yes
Basket	Yes
Wheel	3X3" universal swivel and 2x3" brake
Bearing Capacity	15kg
Max push speed	2m/s
Max tilt	15°
Support arm	Option
Net weight	18.5kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	21kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com