

◆◆ PRODUCT SPECIFICATION ◆◆

MEDIFOCUS


CONFIGURATION	PARAMETER
Model	E08
Column	78*100*810mm
Mounting platform	250*210*5mm
Handle	Yes (210*170*25mm)
Basket	Yes (330*210*135mm)
Humidifier hanger	No
Wire wrap	No
Base	Yes (550*520*50mm)
Wheel	Yes (2X3" universal swivel and 2x3" brake)
Bearing Capacity	15kg
Max push speed	2m/s
Max tilt	15°
Net weight	21.5kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	24kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengtren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

 Email: info@bjmedifocus.com
www.bjmedifocus.com