


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	F02
Column	216*151*452mm(L*W*H)
Mounting platform	
Base	φ600*70mm
Handle	Yes (φ 30*381*517mm)
Basket	
Humidifier hanger	Yes (55*40*16mm)
Wire hanger	
Wheel	Yes (4X4" brake)
Bearing Capacity	50kg
Max push speed	2m/s
Max tilt	15°
Net weight	23kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	25.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tenggren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

