

◆◆ PRODUCT SPECIFICATION ◆◆

MediFocus medical trolley

CONFIGURATION	PARAMETER
Model	F03
Column	210*55*735mm
Base	640*500*37mm
Handle	φ 22*450*160mm
Transfer Basin	800*450*220mm
Humidifier hanger	Yes
Oxygen cylinder holder	150*99*3mm
Strap bracket	120*40*2.5mm
Wheel	4X4" swivel and brake
Bearing Capacity	30Kg
Max push speed	2m/s
Max tilt	15°
Net weight	28kg
Carton dimension	805*635*300mm
Gross weight	30.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.comwww.bjmedifocus.com