


◆◆ PRODUCT SPECIFICATION ◆◆



MediFocus medical trolley for

Lowenstein Leoni ventilator

CONFIGURATION	PARAMETER
Model	G01
Column	Φ 60*893*2.5mm
Mounting platform	295*240*10mm
Handle	170*35*45mm
Basket	NO
Humidifier hanger	55*40*16mm
Base	510*550*30mm
Circuit hanger mounting spot	YES
Wheel	4x4" brake
Bearing Capacity	30Kg
Max push speed	2m/s
Max tilt	15°
Net weight	15.4kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	18kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com