


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	I01
Column	77*77*790mm
Mounting platform	340*285*5mm
Handle	Yes (240*18*6mm)
Basket	Yes (330*210*135mm)
Humidifier hanger	Yes (55*40*16mm)
Oxygen cylinder	Yes 150*99*3mm
Base	Yes (608*535*50mm)
Wheel	Yes (4x4" brake)
Bearing Capacity	30kg
Max push speed	2m/s
Max tilt	15°
Net weight	22kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	26.4kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengen Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

