

◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	K02
Outline size	600*550*910mm
Column	70*135*725mm
Base	600*550*165mm
Basket	330*210*135mm
Humidifier hanger	70*55*15mm
Mounting	350*280*35mm
Handle	Yes
Wheel	4X4" swivel and brake
Bearing Capacity	50kg
Max push speed	2m/s
Max tilt	15°
Net weight	17.4kg
Carton dimension	805*635*300mm (L*W*H)
Gross weight	20kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

 Email: info@bjmedifocus.com
www.bjmedifocus.com
