



◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	K03
Outline	640*500*830mm
Mounting platform	210*172*8mm
Column	210*55*615mm
Base	640*500*37mm
Humidifier hanger	Yes
Wire wrap	Yes
Basket	Yes
Wheel	4X4" swivel and brake
Bearing Capacity	15Kg
Max push speed	2m/s
Max tilt	15°
Net weight	24kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	26.6kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tengen Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com

