


◆◆ PRODUCT SPECIFICATION ◆◆



CONFIGURATION	PARAMETER
Model	M02
Column	Φ 51*610(Φ *H)
Mounting platform	No
Handle	Yes (230*175mm)
Auxiliary bar	Yes (400*25*10mm)
Base	Φ520*40mm
Infusion frame hook	Yes (Φ 10*724*160mm)
Wire hanger	No
Wheel	Yes (3X3" universal swivel and 2x3" brake)
Bearing Capacity	15kg
Max push speed	2m/s
Max tilt	15°
Net weight	9.5kg
Carton dimension	900*570*210mm (L*W*H)
Gross weight	12.1kg

Beijing MediFocus Medical Co.,Ltd

Add: Room 101, building No. 3, No. 11 Tenggren Road, Niulanshan Town, SHUNYI, Beijing, China

Tel: +86 (10) -61796973

Email: info@bjmedifocus.com

www.bjmedifocus.com