

medatro® medical trolley**B12**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: ϕ 600*890mm

Material: Q235 Steel+6063 aluminium

Base size: ϕ 600*70mm

Column size: 78*100*810mm

Platform: 230*245*32mm

IV pole: ϕ 19*780mm

Handle: 230*175mm

Basket: 330*210*135mm

Humidifier hanger: 55*40*16mm

Infusion rod: ϕ 19*780mm

Cable hanger: 204*65*34mm

Castor: 3 inch*5pcs (with 2 brakes)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 10.2kg

Carton dimension: 900*570*210mm

Gross weight: 13kg

