

medatro® medical trolley**B13**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: ϕ 600*990mm

Material: Q235 Steel+6063 aluminium

Base size: ϕ 600*70mm

Column size: 78*100*810mm

Platform: 236*117*10mm

Handle: 230*175mm

Basket: 330*210*135mm

Humidifier hanger: 55*40*16mm

Fixed bar: 400*25*10mm

Castor: 4 inch*5pcs (with 5 brakes)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 11.9kg

Carton dimension: 900*570*210mm

Gross weight: 16kg

