

medatro® medical trolley**B17**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: ϕ 600*890mm

Material: Q235 Steel+6063 aluminium

Base size: ϕ 600*70mm

Column size: 78*100*810mm

Mounting platform: 240*225*8mm

Handle: 170*35*45

Basket: 330*210*135mm

Humidifier hanger: 55*40*16mm

Fixed bar: 300*25*10mm

Castor: 4 inch*5pcs (with 5 brakes)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 12kg

Carton dimension: 900*570*210mm

Gross weight: 15kg

