

medatro® medical trolley**B19**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: ϕ 600*890mm

Material: Q235 Steel+6063 aluminium

Base size: ϕ 600*70mm

Column size: 78*100*810mm

Mounting platform: 100*78*128mm

Handle: 230*175mm

Drawer: 430*290*131.5mm

Basket: 330*210*135mm

Humidifier hanger: 55*40*16mm

Cable hanger: 204*65*34mm

Castor: 4 inch*5pcs (with 5 brakes)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 20kg

Carton dimension: 900*570*210mm

Gross weight: 23kg

