

**medatro® medical trolley****B20**

medical industry mobility solution

**Safety · Durability · Mobility · Customized LOGO**

**Dimension:**  $\phi$ 600\*890mm

**Material:** Q235 Steel+6063 aluminium

**Base size:**  $\phi$ 600\*70mm

**Column size:** 78\*100\*810mm

**Mounting platform:** 300\*230\*8mm

**Handle:** 230\*175mm

**Basket:** 330\*210\*135mm

**Fixed bar:** 400\*25\*10mm

**Humidifier hanger:** 55\*40\*16mm

**Cable hanger:** 204\*65\*34mm

**Castor:** 5 inch\*5pcs (with 2 brakes)

**Load Capacity:** 30kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Net weight:** 14.9kg

**Carton dimension:** 900\*570\*210mm

**Gross weight:** 18kg

