

**medatro® medical
trolley****B23**

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: $\phi 600 \times 1060 \text{mm}$

Base size: $\phi 600 \times 70 \text{mm}$

Column size: $\phi 51 \times 940 \times 3 \text{mm}$

Platform: $\phi 51 \times 55 \text{mm}$

Basket: $330 \times 210 \times 135 \text{mm}$

Handle: $230 \times 175 \text{mm}$

Rob of handle: $\phi 30 \times 222 \text{mm}$

Bracket for column: $80 \times 31 \times 30 \text{mm}$

Castor: 5pcs*3inch (2 with brake)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 10.1kg

Carton dimension: $900 \times 570 \times 210 \text{mm}$

Gross weight: 12.7kg

