

medatro® medical trolley

## B24

medical industry mobility solution

**Safety • Durability • Mobility • Customized LOGO**

**Dimension:** φ600\*890mm

**Base size:** φ600\*70mm

**Column size:** 78\*100\*810mm

**Platform:** 360\*240\*8mm

**Basket:** 330\*210\*135mm

**Humidifier hanger:** 55\*40\*16mm

**Installation plate for humidifier hanger:** 68\*55\*4mm

**Extension plate of circuit hanger:** 500\*35\*18mm

**Installation location of circuit hanger:** φ28\*45mm

**Oxygen cylinder holder:** 150\*99\*3mm

**Bandage bracket:** 120\*40\*2.5mm

**Castor:** 5pcs\*3inch (2 with brake)

**Load Capacity:** 30kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Net weight:** 17kg

**Carton dimension:** 900\*570\*210mm

**Gross weight:** 19.6kg

