

medatro® medical trolley

B26

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: φ600*890mm

Base size: φ600*70mm

Column size: 78*100*810mm

Platform: 215*118*10mm

Basket: 330*210*135mm

Handle: 230*175mm

Rob of handle: φ30*222mm

Humidifier hanger: 55*40*16mm

Installation plate for humidifier hanger: 68*55*4mm

Cable hanger: 204*65*34mm

Castor: 5pcs*3inch (2 with brake)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 8.9kg

Carton dimension: 900*570*210mm

Gross weight: 11.5kg

