

medatro[®] medical trolley

D01

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: φ600*960mm

Base size: φ600*80mm

Column size: φ51*800mm

Platform: 230*245*32mm

Cable hanger: 280*35*2mm

Basket: 330*210*135mm

Handle: 240*18*6mm

Humidifier hanger: 55*40*16mm

Infusion rob: φ19*490*200mm

Castor: 5pcs*3inch (2 with brake)

Load Capacity: 15kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 18.5kg

Carton dimension: 900*570*210mm

Gross weight: 21kg

