

**medatro<sup>®</sup> medical trolley****D04**

medical industry mobility solution

**Safety • Durability • Mobility • Customized LOGO**

**Dimension:** φ600\*960mm

**Base size:** φ600\*80mm

**Column size:** φ51\*800mm

**Humidifier hanger:** 55\*40\*16mm

**Platform:** 400\*314\*5mm

**Basket:** 330\*210\*135mm

**Handle:** 170\*35\*45mm

**Castor:** 5pcs\*3inch (2 with brake)

**Circuit hanger:** Aluminum

**Load Capacity:** 15kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Net weight:** 18.5kg

**Carton dimension:** 900\*570\*210mm

**Gross weight:** 21kg

