

medatro® medical trolley

## E08

medical industry mobility solution

**Safety • Durability • Mobility • Customized LOGO**

**Dimension:** 550\*520\*960mm

**Base size:** 550\*520\*50mm

**Column size:** 78\*100\*810mm

**Platform:** 250\*210\*5mm

**Drawer:** 430\*290\*131.5mm

**Basket:** 330\*210\*135mm

**Handle:** 210\*170\*25mm

**Shelf:** 400\*360\*12mm

**Castor:** 4pcs\*3inch (2 with brake)

**Load Capacity:** 15kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Net weight:** 21.5kg

**Carton dimension:** 900\*570\*210mm

**Gross weight:** 24kg

