

medatro® medical trolley

F02

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: 512*620*600mm

Base size: 512*620*50mm

Column size: 216*151*452mm

Basket: 330*210*135mm

Humidifier hanger: 55*40*16mm

Handle: $\varnothing 30 \times 381 \times 517$ mm

Castor: 4pcs*4inch (with brake)

Support arm: Aluminum

Load Capacity: 50kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 23kg

Carton dimension: 900*570*210mm

Gross weight: 25.6kg

