

medatro® medical trolley

F03

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: 600*550*1070mm

Base size: 640*500*37mm

Column size: 210*55*735mm

Transfer basin: 800*450*220mm

Transfer basin bracket: 920*440*170mm

Oxygen cylinder holder: 150*99*3mm

Bandage bracket: 120*40*2.5mm

Humidifier hanger: 55*40*16mm

Handle: $\phi 22 \times 450 \times 160$ mm

Castor: 4pcs*4inch (with brake)

Support arm: Aluminum

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 28kg

Carton dimension: 805*635*300mm

Gross weight: 30.6kg

