

medatro® medical trolley

H02

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: 581*520*1220mm

Base size: 581*520*50mm

Column size: 1000*122*35.5mm

Platform: 563*225*35mm

Shelf: 500*490*20mm

Handle: \varnothing 25*560*490mm

Monitor stand: \varnothing 30*300*120mm

Castor: 4pcs*4inch (with brake)

Load Capacity: 50kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 37.5kg

Carton dimension: 1300*635*300mm

Gross weight: 42kg

