

medatro® medical trolley**I01**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 608*535*960mm

Base size: 608*535*50mm

Column size: 77*77*790mm

Platform: 340*285*5mm

Handle: 240*18*6mm

Humidifier hanger: 55*40*16mm

Oxygen cylinder holder: 150*99*3mm

Bandage bracket: 120*40*2.5mm

Castor: 4pcs*4inch (with brake)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 805*635*300mm

Net weight: 22kg

Gross weight: 26.4kg

