

**medatro® medical trolley****K03**

medical industry mobility solution

**Safety · Durability · Mobility · Customized LOGO**

**Dimension:** 640\*500\*830mm

**Base size:** 640\*500\*37mm

**Column size:** 210\*55\*615mm

**Platform:** 210\*172\*8mm

**Cable hanger:** 186\*55\*75mm

**Handle:** 240\*18\*6mm

**Humidifier hanger:** 55\*40\*16mm

**Basket:** 330\*210\*135mm

**Castor:** 4pcs\*4inch (with brake)

**Load Capacity:** 15kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Carton dimension:** 900\*570\*210mm

**Net weight:** 24kg

**Gross weight:** 26.6kg

