

medatro® medical trolley**K05**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 600*550*1140mm

Base size: 600*550*165mm

Column size: 70*135*1000mm

Platform: 500*430*30mm

Shelf size: 500*430*30mm

Monitor support: 405*35mm

Drawer: 400*425*95mm

Handle: $\varnothing 22 \times 450 \times 160$ mm

Installing plate for handle: 130*100*8mm

Castor: 4pcs*4inch (with brake)

Load Capacity: 100kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 1060*570*260mm

Net weight: 40kg

Gross weight: 43kg

