

medatro® medical trolley**K06**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 600*550*1140mm
Base size: 600*550*165mm
Column size: 70*135*1000mm
Platform: 500*430*30mm
Shelf size: 500*430*30mm
Monitor support: 405*35mm
Drawer: 400*425*95mm
Handle: $\varnothing 22 \times 450 \times 160$ mm
Installing plate for handle: 130*100*8mm
Castor: 4pcs*4inch (with brake)
Load Capacity: 100kg
Max push speed: 2m/s
Max tilt angle: 15°
Carton dimension: 1060*570*260mm
Net weight: 50kg
Gross weight: 53kg

