

**medatro® medical trolley****K08**

medical industry mobility solution

**Safety · Durability · Mobility · Customized LOGO**

**Dimension:** 600\*550\*980mm

**Base size:** 600\*550\*165mm

**Column size:** 70\*135\*840mm

**Platform:** 500\*430\*30mm

**Shelf size:** 500\*430\*30mm

**Drawer:** 400\*425\*95mm

**Monitor support:** VESA Standard

**Storage box:** 360\*140\*70mm

**Coupling agent pendant:** 130\*65\*50mm

**Cable collecting box:** 780\*100\*77.4mm

**Handle:**  $\varnothing 22$ \*450\*160mm (red/ grey)

**Installing plate for handle:** 130\*100\*8mm

**Castor:** 4pcs\*4inch (with brake)

**Load Capacity:** 100kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Carton dimension:** 670\*620\*970mm

**Net weight:** 49kg

**Gross weight:** 57.6kg

