

medatro® medical trolley**K09**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 600*550*840mm

Base size: 600*550*165mm

Column size: 70*135*725mm

Platform: 500*430*30mm

Shelf size: 500*430*30mm

Monitor support: VESA Standard

Storage box: 360*140*70mm

Coupling agent pendant: 130*65*50mm

Cable collecting box: 660*100*77.4mm

Drawer: 400*425*95mm

Handle: $\varnothing 22$ *450*160mm (red/ grey)

Installing plate for handle: 130*100*8mm

Castor: 4pcs*4inch (with brake)

Load Capacity: 100kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 670*620*970mm

Net weight: 40kg

Gross weight: 47.6kg

