

medatro® medical trolley**K16**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 600*550*980mm

Base size: 600*550*165mm

Column size: 70*135*850mm

Platform: 500*430*30mm

Shelf size: 500*430*30mm

Handle: $\varnothing 22 \times 450 \times 160$ mm (red)

Castor: 4pcs*4inch (with brake)

Load Capacity: 100kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 805*635*300mm

Net weight: 30kg

Gross weight: 33kg

