

**medatro® medical trolley****L01**

medical industry mobility solution

**Safety · Durability · Mobility · Customized LOGO**

**Dimension:** 810\*630\*1160mm

**Base size:** 810\*630\*60mm

**Column size:** 200\*50\*935mm

**Platform:** 610\*500\*30mm

**Shelf size:** 550\*500\*30

**Drawer:** 470\*210\*30mm

**Monitor support:** 60\*330mm

**Handle:** 330\*30\*15mm

**Castor:** 4pcs\*5inch (with 2 brakes)

**Load Capacity:** 100kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Carton dimension:** 950\*760\*1150mm

**Net weight:** 55kg

**Gross weight:** 70kg

